

# U-Step 2

Specifically for  
Neurological  
Conditions

Prevents  
Falls

**MOBILITYCARE**  
*Since 1992*

Phone: (03) 9568 8383  
39 Cambro Road, Clayton VIC 3168  
info@mobilitycare.net.au  
www.mobilitycare.net.au

## The U-Step 2 Walking Stabilizer


.....was designed to increase independence and eliminate falling among those with neurological conditions. The leading-edge features of the U-Step 2, as compared to other walking aids, make it superior by providing more stability, maneuverability and control. Tens of thousands of people with neurological conditions have benefited from the U-Step 2 design.

Regularly prescribed by physicians and physical and occupational therapists from leading hospitals and

rehabilitation centers, and covered by most major insurances and Medicare, the U-Step 2 has been tested and proven to be beneficial for patients with various neurological conditions, including:

- Parkinson Disease • Multiple Sclerosis
- Balance Disorders • Brain Injuries
- ALS • PSP/MSA • Ataxia • Stroke

It's ALL About.....the **BASE!**



Stability...

Our patented U-shaped base, is ultra-stable. It surrounds you and braces you in every direction from falling. You feel as stable as you would feel holding onto another person's arm. Other walkers only partially brace you from falling.

Control...

The advanced control features play an important role in the security that the U-Step provides you. The reversed braking system and the rolling resistance control explained inside are exclusive to the U-Step 2 and put you in complete control.

Maneuverability...

By placing the main turning wheels towards the middle of the base, as opposed to in the rear as in other rolling walkers, the U-Step 2 turns around in 29 inches as compared to 50 inches for most other advanced walkers. This will allow you to turn in tight indoor areas, such as washrooms, kitchens, and hallways. At 23 inches wide, you can easily walk through narrow doorways.



# U-Step<sup>2</sup> Advanced features to address your walking needs

## Reversed Braking System



The innovative braking system is easy to use and puts you in complete control. The **U-Step 2** will not roll until you lightly squeeze either hand brake. Once you release the hand brake, the unit will stop immediately. This feature is particularly helpful when standing up from a chair because the unit will not roll away from you.

## Rolling Resistance Control



Unlike other rolling walkers, you can adjust the rolling speed of the new **U-Step 2**. Many people find that wheeled walkers roll too fast for them. That is why we engineered a control for setting the rolling resistance. Loosen the set screw and rotate the adjustment lever to the desired amount of resistance, then re-tighten the set screw.

## Spring-Loaded Front Wheel



The patented spring-loaded front caster system glides easily over uneven surfaces, such as indoor molding strips and cracks in the sidewalk. This avoids sudden jolts as well as the need to lift the walker over small obstacles.

## Standard Model



Model #: US-PC2

## Laser & Sound Cueing Module



This optional feature is primarily used by those with Parkinson's freezing but can be used by anyone with an irregular gait pattern. This Cueing Module can help get you started, normalize your walking and increase your stride. Press the red button on the module to project a bright red laser line on the floor to guide your step. You can also turn on the sound feature to set a beat pattern for walking speed.

## Easily Transported



Folding the **U-Step 2** is as easy as one, two, three! Lift the release lever by the seat, pull up on the seat support to fold and put it in the car; no bending down to pick it up.

## Comfortable Padded Seat



"I stopped falling! When I bought my **U-Step 2** walker, the results were nothing short of a miracle. I immediately felt safe, as if the walker were a natural extension of my body. It is lightweight, easy to operate and a far more dignified design than other walkers that I tried."  
 — Ralph Bakery, Los Angeles, CA

# 3 Models to meet your specific needs

## Platform Model

Ideal for Stooped Posture, weak upper body, stroke, brain injury. Highly adjustable for optimal setup.



Model #: US-PC2-PL



## Press Down Model

Ideal for those with weak or no hand strength to squeeze a standard hand brake. To go, press down on the brake handle or the flip bar. Choose left or right brake handle.



Model #: US-PC2-PD



# Accessories

## Laser & Sound Cueing Module



Primarily used by those with Parkinson's freezing but can be used by anyone with an irregular gait pattern. The Laser and Sound Cueing Module can help get you started, normalize your walking and increase your stride. Press the red button on the module to project a bright red laser line on the floor to guide your step. You can also turn on the sound feature to set a beat pattern for walking speed.

## Tray & Basket Option



With the tray and basket you can easily carry things at home. Operates on sliding rails, so you can pull it close for easy access and push it away for sitting down. Mounted high for reduced bending, with a large wire basket and tray for carrying plates, drinks, phone, keys and more.



## Why the U-Step was developed

"In 1991, my mother's neurological condition affected her walking to the point of requiring a walking aid. Her physician prescribed a standard walker, in addition to therapy and leg braces. However, the walker hindered her walking, and more importantly did not prevent her from falling. For two years, I worked with physicians, therapists and engineers to develop a better solution. Today, tens of thousands of people with neurological conditions walk more safely and easily with the U-Step 2" says Jonathan Miller, President of In-Step Mobility Products Corporation.

## Quality Assurance

In an effort to create the best product possible, we paid attention to the details of the design and we used only quality components to insure the U-Step 2 would hold up well over time. For example, the ergonomically designed and positioned handle bars reduce stress on the back and joints, and encourage an upright stance. The base of the unit is completely welded, not bolted, so it will not loosen over a period of time.

## Insurance & Medicare Coverage

Insurance companies offer coverage for walkers. The U-Step 2 is an advanced walker which carries a product (HCPC) code of E0147. Medicare generally pays for the U-Step 2 for those with neurological conditions or conditions causing the limited use of one hand.

### U-Step 2 Walking Stabilizer Specifications

Medicare code (HCPC) .....	E0147
Weight capacity .....	375 lbs
Height adjustment range .....	accommodates users 4'10" to 6'2"
Size of Padded Seat .....	19" x 8"
Height from floor to seat (Junior 19") .....	22"
Overall width .....	23"
Turning circle .....	29"
Weight of Standard Model.....	21 lbs
Material of frame .....	Tubular Steel and Aluminum

— Patent Pending —

## In-Step Mobility Products Corp.

8048 Monticello Ave. Skokie, IL 60076

(800) 558-7837 or (847) 676-1275

Fax #: (847) 676-1202

www.UStep.com

e-mail: Walkers@UStep.com