

Oxygen Anytime Anywhere *with the*

InogenOne G3^{HF}



Approved by 45 Airlines
Maintain Oxygen Saturation Level
Small & Compact Design

- ✓ **Approved by 45 Airlines**
- ✓ **45% Smaller, 30% Lighter
27% More Oxygen/Pound
(Compared to INOGEN G2)**
- ✓ **Maintain Correct Oxygen
Saturation Levels**
- ✓ **Designed for 24 Hour Use**

Weight	2.18 kg (Includes Single Battery)
Size	Length: 22cm Width: 7.6cm Height: 21cm
Oxygen Flow	24/7 POC delivery system Five flow settings 1, 2, 3, 4 & 5
Power	AC Power Supply 100-240V, 50-60Hz (auto sensing to allow worldwide use) DC Power Cable: for mobile use in car
Battery	Duration (single battery): Up to 4.7 hours Recharge up to 2-4 hours with AC or DC power Duration (double battery): Up to 10 hours Recharge up to 4-8 hours with AC or DC power
Noise Level	39 dBA (at flow setting 2)
Warranty	3 years limited warranty (Column 1 year warranty)
Operation	Simple control functions
Use	Designed for 24/7 use.

1300 558 947
oxygenolutions.com.au

 **OXYGEN
SOLUTIONS**

InogenOne G3^{HF} System



InogenOne G3^{HF}

The new Inogen One® G3^{HF} delivers the independence of previous Inogen One portable concentrators, but it's smaller, lighter, and makes more oxygen per pound. It's a portable concentrator light enough to be worn care-free by the active oxygen user. That means there are no tanks to refill and you never have to worry about running out of oxygen.

The Inogen One® G3^{HF} is your one solution, your single solution, for oxygen at home or away.



Inogen One®
AC Power Supply*
Item # BA-301



Inogen One® G3^{HF}
DC Power Cable*
Item # BA-306



Inogen One® G3^{HF}
Carry Bag
Item #CA300



Inogen One® G3^{HF}
8 Cell Lithium Ion
Battery*
Up to 4.7 hours run time
Item #BA-300



Inogen One® G3^{HF}
16 Cell Lithium Ion
Battery**
Up to 10 hours run time
Item #BA-316

**Included with system*

***Accessory sold separately*

1300 558 947
oxygenolutions.com.au

 **OXYGEN
SOLUTIONS**