

THINK OF COMFORT

TOPFORMTM

Est. 1965

LUXOR PLATINUM

LUXOR OPTIMA



CHOOSE YOUR SIZE



Petite
(Optima only)



Mini



Medium



Maxi



Tall

FUNCTIONS: Manual Recliner | Electric Recliner | 1 Motor Lift Recliner | 2 Motor Lift Recliner



One Motor



Two Motor

Ultrafresh Antimicrobial Treated Foams



ULTRATHANE MULTI-WAY-STRETCH ON CONTACT SURFACES WHICH IS:

- Extremely Hydrolysis Resistant
- Fluid Proof
- Anti-Bacterial
- Low Water Vapour Transmission
- Fire Retardant
- Hypoallergenic / Non-irritating
- Breathable
- Chemical Resistant for easy cleaning
- Available in Beige or Navy colours

TRIPLE LAYERED PRESSURE MANAGEMENT SEAT FOAM:

- Bottom high density layer provides a firm foundation
- Middle support layer to increase comfort
- Top layer premium memory foam to assist in pressure relief.


LUXURA MEMORY FOAM

- Soft & luxurious feel
- Reduces pressure on the body and improves circulation
- Greater air flow for heat distribution

SIZE	PETITE	MINI	MEDIUM	MAXI	TALL
Seat Height	47cm	47cm	50cm	53cm	54cm
Seat Depth	45cm	49cm	50cm	52cm	53cm
Seat Width	44cm	48cm	48cm	53cm	54cm
Seat To Top Of Back	57cm	60cm	64cm	68cm	74cm
Chair Width	78cm	78cm	78cm	83cm	83cm
Chair Height	99cm	102cm	107cm	117cm	123cm

LUXOR OPTIMA MEDIUM 1M LIFT CHAIR PRESSURE MAPPING RESULTS

UPRIGHT PRESSURE TEST



RECLINED PRESSURE TEST

