

permobil

THE POWER OF MOBILITY™



# JUST THE FACTS...

As the inventors of innovative front-wheel drive technology more than five decades ago, Permobil continues to be the authority with our C300 Corpus® 3G, giving you the stability, agility, and high-caliber performance you expect from us. Its compact size lets you easily navigate tight spaces, while the powerful motors take you anywhere you need to go. The Corpus® 3G Ergo Backrest mimics the shape and comfort of today's Corpus Backrest, but provides more functionality and adjustment. A powerful center-mounted actuator offers up to 50° of tilt, and users can recline up to a full 175°. The legrests are curved for better support, and the armrests feature built-in channels for cable routing, which significantly reduce the need for messy wire ties. Case closed!

**C300**  
CORPUS 3G





### POWER TILT

50° of tilt with two range selections for valuable repositioning.



### LEGREST COMFORT

Adjustable legrests move independently to suit a user's needs.



### POWER RECLINE

Up to 175° of recline to find "just the right level" of comfort.



### INDEPENDENT ARMRESTS

Adjustable armrest height, angle and pitch as well as a convenient flip-up feature.

### AVAILABLE SHROUD COLORS



ONYX BLACK



VOLCANO RED



COBALT BLUE



PURPLE RAIN



PHANTOM GREY



POP STAR PINK

### SPECIFICATIONS

Maximum User Weight	300 lbs	136 kg
Maximum Speed Options	5.5 MPH, 6.5 MPH	8 km/h, 10 km/h
Driving Range <sup>1</sup>	16 miles	25.5 km
Base Width	24"	610 mm
Base Length	40"	1020 mm
Minimum Turning Radius	25"	640 mm
Ground Clearance	3"	75 mm
Optional Seat Elevator	8"	200 mm
Battery Type	Group 34	
Drive Electronics	90A / 120A R-net	
Maximum Slope <sup>2</sup>	10° Up / 6° Down	

Seat to Floor Height	(Elevator)	17.5"	445 mm
	(Seat Tube)	17.5", 18.25", 19"	445, 465, 485 mm
Tilt Options	(Power)	0° – 50°	
	(Manual)	-5° – 45°	
Recline Options	(Power)	85° – 175°	
	(Manual)	85° – 120°	
Legrest Elevation Options	(Power)	85° – 170°	
	(Manual)	85° – 170°	
Armrest Pad Lengths		10", 13", 16"	260, 335, 410 mm
Seat Widths		17", 19", 21", 23"	420, 470, 520, 570 mm
Seat Depths		14" – 22"	370 – 570 mm
Backrest Heights		20", 23" – 28"	500, 570 – 710 mm

<sup>1</sup> ACTUAL DRIVING RANGE WILL VARY BASED ON DRIVING AND BATTERY CONDITIONS.  
<sup>2</sup> REFER TO OWNER MANUAL FOR SAFE DRIVING PRACTICES