

# ENHANCE your Active Lifestyle with added Comfort and Performance



**PATHRIDER 10**

4-Wheel Scooter



**Feather-touch Disassembly**

The sleek sporty **Pathrider® 10** offers added style, comfort and performance for your active lifestyle. Standard features include LED lighting package, backlit battery gauge, feather touch disassembly, wraparound easy-drive tiller, 25.4 cm (10") solid tires and a front basket. A 400 lb. weight capacity, a per charge range of up to 25 km (15.5 miles) and a maximum speed up to 8.85 km/h (5.5 mph) compliment the versatility of the Pathrider 10. Available in Viper Blue.

PRIDE® **MOBILITY SCOOTERS**  
LIVE YOUR BEST®

[pridemobility.com.au](http://pridemobility.com.au)



## Specifications 4W: SC710

Weight Capacity	159 kg (350 lbs.)
Maximum Speed <sup>1</sup>	Up to 8.85 km/h (5.5 mph)
Ground Clearance <sup>2</sup>	7.62 cm (3")
Turning Radius <sup>2</sup>	137.16 cm (54")
Overall Length <sup>2</sup>	119.38 cm (47")
Overall Width <sup>2</sup>	56.5 cm (22.25")
Front Tyres	25.4 cm (10") solid
Rear Tyres	25.4 cm (10") solid
Range Per Charge <sup>2,4</sup>	Up to 25 km (15.5 miles)
Total Weight Without Batteries	62.36 kg (137.5 lbs.)
Heaviest Piece When Disassembled	27.66 kg (61 lbs.) front section
Standard Seating	Type: Foldable, compact, vinyl-covered, moulded plastic Weight: 12.7 kg (28 lbs.) Material: Black Vinyl Dimensions: Width: 46 cm (18") Depth: 43 cm (17")
Motor Output	24-volt DC motor, sealed mini-transaxle
Brakes	Electronic, regenerative and electromechanical
Battery Requirements <sup>3,5</sup>	(2) 12V deep-cycle Size: (2) U-1 (31-36 AH) standard (2) 40 AH optional Weight: (U-1) 11 kg (24.5 lbs.) each (40 AH) 14.74 kg (32.5 lbs.) each
Battery Charger	Off-board, 5 amp
Colour	Viper Blue



## Features

- Exclusive one-hand feather touch disassembly provides the ultimate quick transport system
- Easy-drive tiller with wraparound handles
- Backlit battery gauge
- Exclusive stylish, lightweight, non-scuffing, black, low-profile tires
- Long-lasting, high-intensity LED headlight offers optimal pathway illumination
- Lightweight pinchless hinge seat eight position rotation, sliders and viscoelastic foam inserts for increased comfort
- Easy-grip tiller adjustment knob
- Easy access tiller-mounted charger port
- Auto-connecting front-to-rear harness

## Accessories

- Pride® rear bag (includes crutch/walking stick holder)
- Walker holder (must have Pride® rear bag to attach)
- Oxygen tank holder (must have Pride® rear bag to attach)
- Sun canopy (must have Pride® rear bag to attach)
- Saddlebag
- Mobile phone holder
- Cup holder
- Flag
- Rearview mirror

<sup>1</sup> Varies with user weight, terrain type, battery amp-hour (AH), battery charge, battery condition and tyre condition. These specifications can be subject to a variance of (+/- 10%).  
<sup>2</sup> Due to manufacturing tolerances and continual product improvement, this specification can be subject to a variance of (+ or - 3%).  
<sup>3</sup> AGM or gel cell type required.  
<sup>4</sup> Tested in accordance with ANSI/RESNA, WC Vol2, section 4 & ISO 7176-4 standards. Results derived from theoretical calculations based on battery specifications and drive system performance. Test conducted at maximum weight capacity.  
<sup>5</sup> Battery weight may vary based on manufacturer.

### NOTE: All specifications subject to change without notice.

The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice. Speed and range vary with user weight, terrain type, battery charge, battery condition and tyre pressure.

Theoretical range is calculated under test conditions in accordance with European standards. Actual range in normal use will depend on many factors, including the condition of the vehicle and its batteries, the weight of the driver, correct tyre pressure, ambient temperature, and the gradient and surface of the road or pavement.



20-24 Apollo Drive Hallam, Victoria 3803  
ACN 088 609 661 ABN 66 088 609 661



(Tel) 03 8770 9600 • (Fax) 03 9703 2720 • [www.pridemobility.com.au](http://www.pridemobility.com.au)

**Pride®**  
Mobility Products  
Australia Pty. Ltd.