

# PREMIUM

Dual Motor

*"Premium leather  
recliner"*



CE

The Premium Riser Recliner features a sleek and modern design which will suit any décor and is upholstered in luxurious leather. Its dual-motor mechanism will allow you to operate and move the legrest and the backrest independently.

**drive**<sup>TM</sup>  
[www.drive-medical.com.au](http://www.drive-medical.com.au)

The Premium Riser has been designed with supreme comfort and ease of use in mind, and its attractive Charcoal upholstery makes this chic and stylish chair ideal for any surroundings.



## Superb Comfort, Style & Performance...

- » Dual-Motor chair features two motors, which enable the legrest to be elevated / lowered and backrest to be reclined independently to give more flexible support options.
  - » Modern and stylish backrest design.
  - » It will rise and tilt users to their feet.
  - » Excellent snooze position for supreme comfort and relaxation.
  - » Quiet motors provide a smooth and discreet operation.
  - » Ergonomically contoured seat and back for added comfort and support.
  - » Fully padded chaise for increased comfort.
  - » Quality assured, meeting all relevant fire retardancy standards.
  - » Built-in fail-safe battery back-up system, which enables the chair to be raised and lowered in the event of a power failure\*.
  - » Comes with a simple-to-operate handset.
  - » Stitched handset holder to prevent loss or damage to the handset.
  - » Ideal for individuals who have difficulty standing up and sitting down due to restricted mobility (e.g. those with MS, Arthritis, Rheumatism and circulation problems).
  - » Available in 1 colour option: Charcoal (Leather).
  - » 2-year warranty on wooden frame and electrical components.
  - » Maximum weight capacity 130kg (20st/285lb).
- \* Batteries included.



## PREMIUM

CODE: 884SCHAU

Max Weight Capacity:  
130kg (20st/285lb)



Height:  
116cm (46")



Width:  
76cm (30")



Depth:  
80cm (31")



Seat Depth:  
52cm (21")



Seat Height:  
50cm (20")



Seat Width:  
48cm (18")



### Colour Options



Charcoal

Leaflet Code: LL012

If you would like to view this product please contact your local stockist, details below.