

STELLA

Dual Motor

“Available in micro suede fabric or leather”



CE

The Stella Riser Recliner is an extremely comfortable chair due to its fully padded upholstery, which has been designed in two sections to provide comfortable head and lumbar support.

driveTM
www.drive-medical.com.au

The Stella Riser features a dual-motor mechanism which will allow you to operate and move the legrest and the backrest independently. It is available in either leather or micro-suede fabric, with a choice of 4 colour options.



Superb Comfort, Style & Performance...

- » Two-section padded backrest offering lumbar support and a comfortable headrest.
- » Dual-Motor chair features two motors, which enable the legrest to be elevated / lowered and backrest to be reclined independently to give more flexible support options.
- » It will rise and tilt users to their feet.
- » Excellent snooze position for supreme comfort and relaxation.
- » Quiet motors provide a smooth and discreet operation.
- » Ergonomically contoured seat and back and padded upholstery for added comfort and support.
- » Fully padded chaise for increased comfort.
- » Two convenient side pockets for easy storage.
- » Protective, washable head and armrest covers.
- » Quality assured, meeting all relevant fire retardancy standards.
- » Built-in fail-safe battery back-up system, which enables the chair to be raised and lowered in the event of a power failure*.
- » Comes with a simple-to-operate handset.
- » Stitched handset holder to prevent loss or damage to the handset.
- » Swivel castors enable easy manoeuvrability of the chair.
- » Ideal for individuals who have difficulty standing up and sitting down due to restricted mobility (e.g. those with MS, Arthritis, Rheumatism and circulation problems).
- » Available in 4 colour options: Chocolate and Navy (Leather) or Forest and Mudbrick (Micro-Suede fabric).
- » 2-year warranty on wooden frame and electrical components.
- » Maximum weight capacity 130kg (20st/285lb).

* Batteries included.

STELLA

CODE: 806SCAU (Chocolate Leather)

806NAU (Navy Leather)

806FAU (Forest Micro Suede)

806MAU (Mudbrick Micro Suede)

Max Weight Capacity:
130kg (20st/285lb)



Height:
120cm (47")



Width:
95cm (37")



Depth:
90cm (35")



Seat Depth:
50cm (20")




Seat Height:
47cm (19")




Seat Width:
53cm (21")


Colour Options



Chocolate



Navy



Forest



Mudbrick

Leaflet Code: LLO09 If you would like to view this product please contact your local stockist, details below.

