

Shoprider

889D

ASK FOR THE BRAND SHOPRIDER



SHOPRIDER[®]
mobility

www.shoprider.com.au

RELIABLE MOBILITY GUARANTEED

Shoprider

889D

SHOPRIDER®
www.shoprider.com.au *mobility*

The SHOPRIDER™ 889D is scooter designed to assist operators who need a larger width scooter than is normally available. The scooter is well featured, it offers exceptional seating capabilities that includes a sliding seat for added flexibility. The 889 D has great features as standard including front and rear suspension, adjustable steering tiller, rear indicators, front light, sliding bench seat and fold up arm rests. The reliability of this impressive scooter by SHOPRIDER™ is further enhanced by the rear wheel drive and electromagnetic brake.

SPECIFICATIONS

889D

Length	1480mm
Width	910mm
Height	980mm
*Safe Climbing Angle	9 Degrees
*Max Climbing Angle	12 Degrees
*Maximum Load <u>Single user</u>	225Kgs
Charger	8Amp
Motor Size	2.0Hp
Maximum Speed	12Kph
Battery capacity min	2 x 75 Amp
Maximum range#	up to 26km
Wheel Size	330mm
Turning Radius	2550 mm
Weight (with battery)	159kg

FEATURES:


- Wide bench Seat
- Loop Steering Tiller
- Adjustable Steering Tiller
- Adjustable width arm rests
- Sliding Seat
- Front / Rear suspension
- Battery Power indicator
- Head Lights
- Rear Lights & Indicators
- Hazard Lights
- Shopping Basket
- Comprehensive Warranty**

WARNING: Do not exceed 10kph on public roads or paths

IMPORTANT: Always be sure to read your owner's manual thoroughly before operating your machine.

Driving range varies with the factors shown below

- (1) The weight of the occupant
- (2) Level ground & Ambient temperature
- (3) Battery condition / Type of charger
- (4) Driving Style



*(Single user)

Shoprider is listed and complies with the regulations of the Therapeutic Goods Association of Australia as well as the regulatory bodies and standards in America (FDA) and Europe (CE)

SHOPRIDER®
www.shoprider.com.au *mobility*

Your Authorised dealer

For your nearest Distributor call 1800 672 633
Or visit our Web Site at www.shoprider.com.au