





The Scout is ideal for use indoors and outdoors where the terrain is forgiving (eg. smooth paths, shopping centres, etc). The Scout has been developed to ensure a focus on Comfort, Style, Performance and Reliability.



Superb Comfort, Style & Performance...

- Top speed of 4mph*
- Maximum range of 10 miles (12ah) and 13 miles (20ah) on full battery charge*
- Comfortable lightweight seat which swivels 360°
- Removable battery pack with easy to carry handle and locking mechanism
- On board and in-house charging as standard
- Automatic electromagnetic brake system
- > Tool free height adjustable seat
- > Delta Bars as standard
- Easy to change body panels

- > Easy tiller adjustment
- Removable seat, fold down tiller and splits easily for storage and transportation
- > Stylish silver wheel hubs
- Black puncture proof non-marking tyres
- > Width adjustable armrests
- Front removable storage basket as standard
- > Rear anti tip wheels as standard
- > 12ah & 20ah batteries
- Weight capacity 130kg (20st)
- * Depends on user mass, terrain and battery condition.



Splits easily for storage and transportation.

CODE: MS009 (4 - Wheels 20 ah) CODE: MS008 (4- Wheels 12 ah) CODE: MS003 (3- Wheels 12 ah)

DESCRIPTION: Scout



Length: 108cm (42.5")



Height: 91cm (36")



Width: 48cm (19")



Max Capacity: 130kg (20st)



Colour Options

Weight: 43kg (95lbs)



Max Gradient:



Delta Bars as standard.



Splits easily in seconds.





Leaflet Code: LL170

If you would like to view this product please contact your local stockist, details below.



360° swivelling seat.



Removable battery.



1300 008 267

http://ilsau.com.au/product/scout-mobility-scooter/