

# BELLAVITA

BathLift

## bellavita

*"A true bathing experience"*



CE

This ultra modern bathlift takes the Bellavita away from an austere and clinical look and is more in keeping with today's requirements. Developed with the help of Nursing Professionals, it includes several excellent new features. Reaches a record seat height of 48cm. The lightest bathlift on the market at only 9.3kg. Backrest reclines to a market leading 50°. At 6cm the Bellavita is the lowest bathlift available.

**drive**<sup>TM</sup>  
[www.drivemedical.co.uk](http://www.drivemedical.co.uk)

This ultra modern bathlift takes the Bellavita away from an austere and clinical look and is more in keeping with today's requirements. Developed with the help of Nursing Professionals, it includes several excellent new features...



Lightweight, fully waterproof floating hand control



The lightest bathlift on the market at only 9.3kg



## For Those Who Require Their Bathlift To Be Feature Rich...

- » The lightest bathlift on the market at only 9.3kg
- » At 6cm the Bellavita is the lowest bathlift available
- » Backrest reclines to a market leading 50 degrees
- » Reaches a record seat height of 48cm
- » Lifetime guarantee on frame & motor
- » Ergonomically designed for greater comfort and safety to the user
- » Easy to install and operate
- » User friendly charger
- » Quick charge lithium battery technology
- » Built in safety controls
- » Lightweight, fully waterproof floating hand control
- » Washable hygienic covers available in white, blue or grey
- » Hygienic cutaway and slimline profile for easy handling
- » Easy to dismantle for better storage and transportation
- » 4 quick release suckers with option for an additional 2 for anti-slip bath base
- » High weight capacity - 22 stone (140kg)
- » 2 years guarantee on battery/handset
- » Headrest available as an optional extra

## TECHNICAL DATA

Overall length (with backrest in a reclined position):	830 mm
Overall width (seat incl. side flaps):	700 mm
Width of seat:	350 mm
Depth of seat:	420 mm
Overall length of seat:	570 mm
Sitting height min:	60 mm
Sitting height max:	480 mm
Height of backrest:	640 mm
Width of backrest:	320 - 360 mm
Backrest reclines:	50°
Width of baseplate:	310 mm
Length of baseplate:	560 mm
Load capacity max:	140kg
Total weight (excl. hand control):	9,3kg (20.5 lb)

Leaflet Code: LL127

If you would like to view this product please contact your local stockist, details below.

