

IC333

Hi-Lo Trendelenburg Adjustable Bed



Head Lift Angle	62°
Foot Lift Angle	42°
Trendelenburg Angle	15°



HEIGHT ADJUSTABLE



TRENDELENBERG TILT



HEAD & KNEE



200KG WEIGHT LIMIT




10 YEAR WARRANTY

- Dual lift motors totalling 12,000 Newtons lifting power
- Heavy duty locking castors
- Lowest point to top of base when lowered is 22cm
- Highest point to top of base when raised is 66cm
- Available in Onyx & Stone fabric




Long Single

Base:
(2080 x 960)
Mattress:
(2030 x 900)




King Single

Base:
(2080 x 1120)
Mattress:
(2030 x 1070)



Double 1 Piece

Base:
(2080 x 1400)
Mattress:
(2030 x 1350)



Queen 1 Piece

Base:
(2080 x 1580)
Mattress:
(2030 x 1500)



King (per side)

Base:
(2080 x 960)
Mattress:
(2030 x 900)



SLEEPSYSTEMSAUSTRALIA
INNOVATORS IN SLEEP

i-care
Caring for your comfort

IC333 Headboards

- Upholstered Headboard to match IC333
- Available in all bed sizes
- Easily attachable and detachable
- 5 year warranty



Companion Base

Non Adjustable

- Non mechanical upholstered slat base for caregivers
- Available in onyx and stone fabric



Long Single




King Single

Base:
(2080 x 960)
Mattress:
(2030 x 900)

Base:
(2080 x 1120)
Mattress:
(2030 x 1070)

Choice of Mattress



IC20 - Visco Medium

- 20cm high density visco elastic and support foam
- Suitable for low to medium risk pressure care



IC25 - Visco Plush

- 22.5cm high density visco elastic and support foam
- Suitable for low to medium risk pressure care



Euroflex Latex - Medium

- 18cm of natural ventilated latex
- Suitable for low risk pressure care

Accessories




Bedstick

- Easy removable from brackets
- Height: 91cm



High Side Rail

- Easy removable from brackets
- Height: 90cm Width: 36cm



Low Side Rail

- Easy removable from brackets
- Height: 58cm Width: 68cm



Full Length Fold-down bed Rail

- Fold down mechanism
- Height: 44cm Width: 138cm



Over Bed Table

- Height adjustable
- Tilt adjustable



Over Bed Pole

- Free standing
- Adjustable handle height



SLEEPSYSTEMSAUSTRALIA
INNOVATORS IN SLEEP

i-care
Caring for your comfort