

## WALKBELT

Slip resistant and soft padded

- ✓ Assists carer to safely transfer patient in comfort.
- ✓ Heavy duty fabric for maximum safety.
- ✓ Double stitched, highly visible hand grips.
- ✓ Slip resistant inner lining to minimise belt slipping.
- ✓ Improved edging to provide longevity.



**INSTRUCTION:**

- Select correct size
- Ensure buckle is securely fastened.
- Hand wash with a mild detergent only.
- Hang dry.

**SMALL/MEDIUM:**

Minimum length 69cm  
Maximum length 114cm

**LARGE/XLARGE:**

Minimum length 89cm  
Maximum length 144cm

**CAUTION:**

- Monitor weakness or malfunction in the buckle.  
(If caused by excessive use)
- Ensure buckle is securely connected.
- Check for loose stitching on the handle grips.