EASY RIDER

Scooter



Inspired by the success of the Sport Rider, the design team got to work on a new model that would continue to redefine the typical luxury scooter.



The Easy Rider has a striking new design, stunning instrumentation and offers precise handling that must be experienced to be appreciated.




drive

Telescopic front motorcycle suspension for optimal handling



Motorcycle style wheels with chrome mud guards



Twin dial digital dash



All round lighting

Superb Comfort, Style & Performance...

- » Top speed of 8mph*.
- » Maximum range of 31 miles on full battery charge*.
- » Comfortable rotating sport style seat with sliding and reclining facility.
- » Height adjustable headrest and angle adjustable armrests to suit user requirements.
- » Angle adjustable handle bars allows comfortable positioning for the user (requires tools).
- » Throttle lever with automatic safe-braking system.
- » Additional Hand Brake for added security.
- » Telescopic front motorbike suspension for optimal handling.
- » Illuminated controls feature a twin dial digital dash.
- » Trip computer with odometer for measuring distance travelled.
- » Easy to use controls.
- » All round lighting, including large headlight, front/rear indicators and rear brake lights to improve outdoor visibility at dusk or when dark.

- » Twin adjustable rear view mirrors for extra visibility.
- » Large accommodating footwell with anti-slip rubber mat.
- » Motorbike style wheels with modern chrome mud guards.
- » Rear bumper provides additional protection.
- » Rear anti tip wheels as standard.
- » Lap belt included for added safety.
- » Large seat back storage pocket for keeping essentials safe when out and about.
- » Anti roll back system safety device prevents scooter from rolling backwards on hills.
- Freewheel facility allows movement of scooter without switching the motor on.
- » Optional lockable and removable rear pannier box.
- » Recommended MK batteries: 2 x 12V 79Ah or 2 x 12V 96Ah.
- » Weight capacity 180kg (28st).
- * Subject to conditions of use.

CODE: SROO3BLK DESCRIPTION: Easy Rider



Length: 166cm (65")



Height: 115cm (45")



Width: 76cm (30")



Max Capacity: 180kg (28st)



Ground Clearance: 12cm (5")



Max Gradient:

Leaflet Code: LL209

If you would like to view this product please contact your local stockist, details below.