

EASY RIDER

Scooter



“State of the art
performance and design”



The Easy Rider has a striking new design,
stunning instrumentation and offers
precise handling that must be
experienced to be appreciated.

www.drive-medical.com.au

drive
DeVilbiss
HEALTHCARE



EASY RIDER

Scooter

Superb Comfort, Style & Performance...

- Top speed up to 15kph (9mph)*.
 - Maximum range up to 50km (31 miles) on full battery charge*.
 - Comfortable rotating sport style seat with sliding and reclining facility.
 - Height adjustable headrest and angle adjustable armrests to suit user requirements.
 - Angle adjustable handle bars allows comfortable positioning for the user (requires tools).
 - Throttle lever with automatic safe-braking system.
 - Additional Hand Brake for added security.
 - Telescopic front motorbike suspension for optimal handling.
 - Illuminated controls feature a twin dial digital dash.
 - Trip computer with odometer for measuring distance travelled.
 - Easy to use controls.
 - All round lighting, including large headlight, front/rear indicators and rear brake lights to improve outdoor visibility at dusk or when dark.
 - Twin adjustable rear view mirrors for extra visibility.
 - Large accommodating footwell with anti-slip rubber mat.
 - Motorbike style wheels with modern chrome mud guards.
 - Rear bumper provides additional protection.
 - Rear anti tip wheels as standard.
 - Lap belt included for added safety.
 - Large seat back storage pocket for keeping essentials safe when out and about.
 - Anti-roll back system safety device prevents scooter from rolling backwards on hills.
 - Freewheel facility allows movement of scooter without switching the motor on.
 - Optional lockable and removable rear pannier box.
 - Recommended MK Gel batteries: 2 x 12V (97Ah)**.
 - Weight capacity 180kg (28st/392lb).
- * Subject to conditions of use.
** Batteries not included



Telescopic front motorcycle suspension for optimal handling



Motorcycle style wheels with chrome mud guards



Twin dial digital dash



All round lighting

CODE: SR003BLKAU
DESCRIPTION: Easy Rider



Length:
166cm (65")




Height:
115cm (45")



Width:
76cm (30")



Max Capacity:
180kg (28st/392lb)



Ground Clearance:
12cm (5")



Max Gradient:
10°