



SKU: S11100

ABLE LIFE UNIVERSAL STAND ASSIST

\$299.99

The able life universal stand is a perfect for any couch, reclinable lift chair or arm chair. The stand-assist provides dual support with its cushion handles that aids you to get to a standing position. It has a universal fit and can be fitted on top the required chair or in between chair cushions when necessary. The stand-assist is fully adjustable in both height and depth.



Scan QR code to
view product online

PRODUCT DESCRIPTION

Features & Benefits:

- Dual Support Handle - makes sitting and standing easy
- Adjustable -adjusts in height to fit most couches/chairs
- Easy installation - snap together, slide it into the couch and adjust the height
- Low profile - matches look and feel of your home
- Powder coated steel construction

ADDITIONAL INFORMATION

Brand	Able Life
Handle height range from floor	54 – 80 cm
Floor to base height	35 – 78 cm
Handle height adjustment from base	19 – 26 cm
Depth	60 – 68 cm
Width adjustment	48 – 65 cm
Weight of product	3 kg
Weight capacity	136 kg