



SKU:

INVACARE ACTION 3NG LEVER DRIVE WHEELCHAIR

\$3,999.99

The Invacare Action 3NG wheelchair has been re-engineered to make sure you receive a superior driving performance. The Invacare 3NG is designed with a lightweight frame and high quality mechanisms that makes it much easier to roll and maneuver the wheelchair. The wheelchair's unique one arm lever drive control system gives you the ability to propel the chair easily. Using just one arm through the lever aids to propel the chair forward or backward, and its ability to steer and brake has been integrated carefully to ensure your safety.

The Invacare Action 3NG'S Benefits include,

- Versatile adjustment options
- Action3 offers rapid and easy assessment
- The swing in and out leg rest is equipped with an easily operated quick release mechanism
- Elevating leg rests offers optimal length compensation
- With a total weight of 14kg, the 3NG is easy to dismantle for transport in a small car

- Enables easy handling and ensures improved stability
- With its modular design the 3NG is easy to refurbish with the range of options provided
- The footrest swings in and out for easy transfers and has a quick release mechanism
- Rear wheel positions can be adjusted and set to Standard, Medium, Active or Amputee positions



Scan QR code to
view product online

VARIATIONS








Image	SKU	Price	Description	Sizes	Type
	S8737	\$3,999.99		16 x 16	Left
	S12553	\$3,999.99		18 x 16	Left
	S2731	\$3,999.99		18 x 18	Left
	S10756	\$3,999.99		20 x 18	Left

Image	SKU	Price	Description	Sizes	Type
	S12552	\$3,999.99		16 x 16	Right
	S9550	\$3,999.99		18 x 18	Right
	S10579	\$3,999.99		20 x 18	Right

PRODUCT DESCRIPTION

The Invacare Action 3NG has been re-engineered to make sure users receive a superior driving performance. The 3NG wheelchair is equipped with a lightweight frame and high quality mechanisms that makes it much easier for users to roll and maneuver the wheelchair. The wheelchair's one arm lever drive control system gives the user ability to propel the chair easily. Using just one arm through the lever aids to propel the chair forward or backward, and its ability to steer and brake has been integrated carefully to ensure safety.

ADDITIONAL INFORMATION

Brand	Invacare
Max SWL	125kg (275 lbs)
Total width	55cm - 68cm
Seat width	40cm (16") / 45cm (18")
Total height	80 - 100cm
Total length	98cm
Seat height	42.5cm
Item weight	14.5kg (31.9 lbs)
Sizes	16 x 16, 18 x 16, 18 x 18, 20 x 18