

**SKU:** S10425




## **AIRSEP FOCUS REPLACEMENT BATTERY**

**\$349.99**

The newly designed Airsep Freestyle has an easy access front to swap battery packs on the go. The Airsep Freestyle has a re-designed Power Supply Connection allowing multi-directional installation with no pins to align and no release button to press making it simple and easy to use.

The battery pack will re-charge anytime AC or DC power is connected. Charge time is approximately 4 hrs.



Scan QR code to  
view product online

## ADDITIONAL INFORMATION

**Brand**

AirSEP