

SKU: S10425




AIRSEP FOCUS REPLACEMENT BATTERY

\$349.99

The newly designed Airsep Freestyle has an easy access front to swap battery packs on the go. The Airsep Freestyle has a re-designed Power Supply Connection allowing multi-directional installation with no pins to align and no release button to press making it simple and easy to use.

The battery pack will re-charge anytime AC or DC power is connected. Charge time is approximately 4 hrs.



Scan QR code to
view product online

ADDITIONAL INFORMATION

Brand

AirSEP