




SKU:

ALL PURPOSE SLING

The All Purpose Sling is a 2-in-1 sling designed to give greater flexibility for client use. Manufactured by Tech-Assist, it incorporates head support to allow for a full body lift but can also be used without the head support as a general purpose style sling by folding the head section out of the way. The head section foam can also be removed for comfort if needed. Webbing attachments on the sling are covered to prevent chafing. Fully washable and quick drying.



Scan QR code to
view product online

PRODUCT DESCRIPTION

All Tech-Assist slings are available for children and adults and in many sizes, weight classes and materials, including Polyester Mesh, Solid Polyester, Nylon, or Spacer Mesh (mesh and a wetsuit-like material all-in-one).

Tech-Assist slings have been tested to the requirements of ISO/ARTG: AS ISO 10535 2011.

ADDITIONAL INFORMATION

Maximum weight limit	250 kg
-----------------------------	--------