




SKU: S13166

BARIATRIC OVER TOILET AID EXTRA WIDE - ALUMINIUM RUST FREE

\$199.99

The Hero Medical over toilet aid comes with a robust aluminium lightweight frame (4kg) with a 53cm width between the arms. The seat and lid are made from a durable, high-quality moulded plastic with minimal flexibility. This product has been treated for corrosion resistance.



Scan QR code to
view product online

PRODUCT DESCRIPTION

The Hero Medical over toilet aid has height adjustable legs allowing users to alter the seat height and safely transfer on and off the toilet seat. Raised armrests provide additional support for users when sitting and standing. This toileting aid has a SWL of 205kg.

- Height adjustable legs
- Raised armrests for additional support
- Robust seat and lid
- Durable anti-slip rubber tips, providing excellent stability
- SWL 205Kg
- 53cm width between the arms
- Frame has been treated with corrosion resistance and powder coated paint

ADDITIONAL INFORMATION

| | |
|---|--------------|
| Brand | Hero Medical |
| Material | Aluminium |
| Overall Width | 63cm |
| Overall Depth | 55cm |
| Seat Width | 34cm |
| Width Between the Arms | 53cm |
| Seat Height Range | 41cm – 54cm |
| Armrest Height (from seat to top of arm) | 22cm |
| Unit Weight | 4kg |
| Weight Capacity / Max SWL | 205kg |