



**SKU:**

## **BED PAD WITH TUCK INS**





**\$0.00**

This Bed Pad with Tuck-Ins and a breathable, waterproof polyurethane backing has great capacity to absorb up to 2.5 litres of body fluids. The Bed Pad with Tuck Ins consists of two layers of fabric. Where in contact with the skin it is made of brushed polyester which has been quilted to a lower layer of reinforced needled rayon. Light Blue in colour.



Scan QR code to  
view product online

VARIATIONS

Image	SKU	Price	Description	Size
	S206	\$69.99		Single
	S12045	\$89.99		Double
	S12046	\$74.99		King Single
	S12044	\$99.99		Queen

## **PRODUCT DESCRIPTION**

## ADDITIONAL INFORMATION

**Width**

Single 900 mm  
Double 900 mm  
Queen 900 mm  
King Single 900 mm

**Length**

Single 1000 mm  
Double 1420 mm  
Queen 1540 mm  
King Single 1100 mm

**Absorption Capacity**

2.5 Litres

**Manufacturer's part number**

3058