



SKU:


BED RAISER

\$0.00

This bed raiser is a tough plastic brick for raising beds and chairs, for easier access. Raises furniture by 40mm and 100mm. Sold as each.

Height furniture is raised per unit: 40 mm or 100mm Diameter of well in top: 65 mm Depth of well in top: 15 mm Allowable weight, set of four: 250 kg Material: Polypropylene Cleaning methods: Autoclavable, or use common soaps, detergents

Rubber mats sold separately



Scan QR code to
view product online

VARIATIONS

Image	SKU	Price	Description	Size
	S1215	\$10.99		40mm
	S7269	\$14.29		100mm

ADDITIONAL INFORMATION

Height furniture is raised per unit	40 mm or 100mm
Diameter of well in top	65 mm
Depth of well in top	15 mm
Allowable weight	250 kg
Material	Polypropylene