




SKU: S2173

BEDSIDE COMMODE

\$249.99

This Bedside Commode can be stacked up to 10 high and best suits a care environment. It has a clip on toilet seat, which clips on in a matter of seconds. There are no metal fittings on the toilet seat and this commode is height adjustable and easy to clean. It has a light blue deluxe padded overseat, as well as a removable round bowl with lid and handle.



Scan QR code to
view product online

PRODUCT DESCRIPTION

The commode's lid helps to eliminates spillages and reduce the risk of cross infection. The bedside commode is finished in a powder coated dark blue frame to combat rust and corrosion.

ADDITIONAL INFORMATION

Max SWL	125 kg
Height	830 mm – 980 mm
Width	550 mm
Depth	550 mm
Seat Height (padded)	485 mm – 635 mm
Seat Height (aperture)	450 mm – 600 mm
Seat Width	450 mm
Seat Depth	430 mm
Floor to Top of Armrest	640 mm
daily-living	Commodes
Brand	Days Healthcare