



SKU: S11499

BELLAVITA BATH LIFT

\$999.99

Exclusive to ILS, the Bellavita Bath Lift enables safe and comfortable bathing for people with limited mobility. Its lightweight, battery-powered design offers easy push-button operation and gentle recline, with smooth lowering and raising for confident transfers.



Scan QR code to
view product online

GALLERY IMAGES



PRODUCT DESCRIPTION

The Bellavita Bath Lift is a reliable, lightweight solution that allows individuals to bathe safely and independently. With easy-touch controls, it lowers the seat to just 6 cm, enhancing immersion and comfort. After bathing, it gently lifts the person back up for a safer transfer.

This ergonomic lift reclines to 50°, supporting those who need a wider hip-to-back angle or have restricted flexibility. At a raised height of 48 cm, it aligns with standard wheelchair or commode heights, allowing for level transfers.

The unit splits into two lightweight parts for convenient handling and transport.

A waterproof, floatable hand control features large buttons and built-in safety indicators. The lift will not lower unless there is enough charge to raise it again, while an emergency stop button provides added reassurance.

Four quick-release suction feet secure the lift in the bath, and washable covers support hygiene and easy maintenance.

Who it's for

- For clients with limited lower limb strength, balance challenges or restricted hip movement which impacts safer transfers in and out of the bath.
- For therapists prescribing assistive technology to enhance safe, independent bathing at home or in community care settings.

Key features

- **Push-button lift and lower function** – supports independent, strain-free transfers
 - **50° backrest recline** – ideal for those needing a wider hip-to-back angle
 - **Ultra-low 6 cm seat height** – promotes full immersion, even in shallow baths
- **Raised seat height of 48 cm** – aligns with wheelchairs or commodes for level transfers
 - **Two-part, lightweight design** – easy to transport, install and store
- **Four quick-release suction feet** – ensure secure placement and easy removal
- **Floatable, waterproof hand control** – straightforward operation in or out of the water
- **Lithium battery with charge indicator** – prevents lowering if insufficient power to raise
 - **Emergency stop button** – adds safety and peace of mind
- **Durable, recyclable plastic construction** – With a safe working load up to 140 kg

ADDITIONAL INFORMATION

Safe Working Load (SWL)	140 kg
Gross Weight	9.3 kg
Overall Width	70 cm
Overall Length	83 cm
Seat Width	35 cm
Seat Depth	42 cm
Overall Seat Length	57 cm
Seat Height	6 cm – 48 cm
Backrest Height	64 cm
Backrest Width	32 cm – 36 cm
Backrest Reclines	50°
Baseplate Width	31 cm
Baseplate Length	56 cm
Brand	Drive Medical

[Drive Medical Bellavita Bath Lift User Manual](#)

[Drive Medical Bellavita Bath Lift Brochure](#)

Talk to us on the phone

Our experienced phone sales team can assist you with any product enquiries and help you make a purchase directly on the phone. We are available anytime between 9am and 5pm, Monday to Friday.

[1300 008 267](tel:1300008267)

HOW HEAVY IS THE BELLAVITA BATH SEAT?

The Bella Vita Auto Bath Seat is the lightest bathlift on the market at only 9.3kg.