



SKU: S11499

BELLAVITA BATH LIFT

\$999.99

EXCLUSIVE TO ILS.

Enjoy a bath safely and comfortably with the #1 rated bath lift in UK, the Drive Medical Bellavita reclining bath lift.

Made in Germany, the Bellavita bath lift is the lightest battery-powered reclining bath lift available, weighing only 9.3kg.

Developed with the help of nursing professionals, it is ergonomically designed for you to enjoy a great soak in the bath with a record seat height of 48cm, very low seat travel of 6cm and a backrest that reclines to a market leading 50 degrees.

Built with many safety features, this ultra-modern bath lift is much loved by our customers as a daily living aid.



Scan QR code to
view product online

PRODUCT DESCRIPTION

Enjoy a bath with confidence using the lightest battery powered bath lift available - the Bellavita bath lift. At the touch of a button using the hand control of this ergonomically-designed bath lift, users with restricted mobility are gently lowered down into the bath. The bath lift also lifts the user to the top again for safely getting out.

Key features include:

- At the push of a button, it lowers persons with limited mobility to the bottom of the tub for bathing and relaxation, then elevates them back to the top of the tub for a safe exit.
- The backrest reclines to a market-leading 50 degrees to accommodate customers requiring a greater hip-to-back sitting position, such as those with stiff hips.
- The cutout in the seat's front facilitates personal hygiene.
- Can be used with or without covers; washable hygienic covers are simply fitted/removed and offer additional bathing comfort (available in white or blue)
- Extremely low seat travel of 6 cm (2.4") from the bottom of the bath enhances water immersion and is perfect for shallow baths.
- Record seat height of 48 cm (18.9") allows for safe level transfers in deeper baths
- Weighs only 9.3 kg (20.5lb) to make transporting and handling easier.
- Quickly and simply separates into two lightweight halves for simple installation/removal, storage, and travel
- Contoured, slim-line profile facilitates handling and provides a contemporary, aesthetically pleasant appearance that complements most bathrooms.
- Four quick-release suction feet permit combat-free removal from bath
- Lightweight, completely waterproof, floatable hand control with easy-to-use buttons
- The handheld remote's lithium-battery technology allows for convenient recharging in numerous locations.
- An indicator light illuminates when charging is necessary, and a safety mechanism stops the bath lift from lowering without sufficient energy to return it to the top of the bath.
- The emergency stop button built into the handle ensures maximum safety
- User-friendly jack plug connection to charger for simplicity of usage
- Produced from high-quality, lightweight plastic that is recyclable and simple to clean
- Extensive array of available accessories to optimise bath installation and serve a variety of consumers
- Weight limit of 140 kg

Additional item can be purchased separately:

- [Washable hygienic covers](#) – Available in blue, the covers are easily fitted/removed and provide extra comfort when bathing.

ADDITIONAL INFORMATION

Brand	Drive Medical
Max user weight (SWL)	140kg
Total weight (excl. hand control)	9.2kg
Overall length	83cm (32.5 inches) – with backrest in a reclined position
Overall width	70cm (27.5 inches) – seat including side flaps
Seat width	35cm (14 inches)
Seat depth	42cm (16.5 inches)
Seat length	57cm (22.5 inches)
Seat height	6cm – 48cm (2.5 – 19 inches)
Backrest height	64cm (25 inches)
Backrest width	32 – 36cm (13 – 14 inches)
Backrest reclines	50 degrees
Baseplate width	31cm (12 inches)
Baseplate length	56cm (22 inches)
Warranty	2 years warranty on battery/handset.
Colour	Blue