



**SKU:** S11499

## **BELLAVITA BATH LIFT**

**\$999.99**

Exclusive to ILS, the Bellavita Bath Lift enables safe and comfortable bathing for people with limited mobility. Its lightweight, battery-powered design offers easy push-button operation and gentle recline, with smooth lowering and raising for confident transfers.



Scan QR code to  
view product online

### **GALLERY IMAGES**





## PRODUCT DESCRIPTION

The Bellavita Bath Lift is a reliable, lightweight solution that allows individuals to bathe safely and independently. With easy-touch controls, it lowers the seat to just 6 cm, enhancing immersion and comfort. After bathing, it gently lifts the person back up for a safer transfer.

This ergonomic lift reclines to 50°, supporting those who need a wider hip-to-back angle or have restricted flexibility. At a raised height of 48 cm, it aligns with standard wheelchair or commode heights, allowing for level transfers.

The unit splits into two lightweight parts for convenient handling and transport.

A waterproof, floatable hand control features large buttons and built-in safety indicators. The lift will not lower unless there is enough charge to raise it again, while an emergency stop button provides added reassurance.

Four quick-release suction feet secure the lift in the bath, and washable covers support hygiene and easy maintenance.

### Who it's for

- For clients with limited lower limb strength, balance challenges or restricted hip movement which impacts safer transfers in and out of the bath.
- For therapists prescribing assistive technology to enhance safe, independent bathing at home or in community care settings.

### Key features

- **Push-button lift and lower function** - supports independent, strain-free transfers
  - **50° backrest recline** - ideal for those needing a wider hip-to-back angle
  - **Ultra-low 6 cm seat height** - promotes full immersion, even in shallow baths
- **Raised seat height of 48 cm** - aligns with wheelchairs or commodes for level transfers

- **Two-part, lightweight design** - easy to transport, install and store
- **Four quick-release suction feet** - ensure secure placement and easy removal
- **Floatable, waterproof hand control** - straightforward operation in or out of the water
- **Lithium battery with charge indicator** - prevents lowering if insufficient power to raise
  - **Emergency stop button** - adds safety and peace of mind
- **Durable, recyclable plastic construction** - With a safe working load up to 140 kg

## ADDITIONAL INFORMATION

<b>Safe Working Load (SWL)</b>	140 kg
<b>Gross Weight</b>	9.3 kg
<b>Overall Width</b>	70 cm
<b>Overall Length</b>	83 cm
<b>Seat Width</b>	35 cm
<b>Seat Depth</b>	42 cm
<b>Overall Seat Length</b>	57 cm
<b>Seat Height</b>	6 cm - 48 cm
<b>Backrest Height</b>	64 cm
<b>Backrest Width</b>	32 cm - 36 cm
<b>Backrest Reclines</b>	50°
<b>Baseplate Width</b>	31 cm
<b>Baseplate Length</b>	56 cm
<b>Brand</b>	Drive Medical

[Drive Medical Bellavita Bath Lift User Manual](#)

[Drive Medical Bellavita Bath Lift Brochure](#)

### **Talk to us on the phone**

Our experienced phone sales team can assist you with any product enquiries and help you make a purchase directly on the phone. We are available anytime between 9am and 5pm, Monday to Friday.

[1300 008 267](tel:1300008267)

### **HOW HEAVY IS THE BELLAVITA BATH SEAT?**

The Bella Vita Auto Bath Seat is the lightest bathlift on the market at only 9.3kg.