




SKU:

SUNRISE MEDICAL BREEZY BASIX 2 FIXED BACK FOLDING WHEELCHAIR





~~\$799.99~~ \$799.99

The Breezy BasiX² is a lightweight folding wheelchair that features amazing adjustability and rigidity, combined with a modern and clean look. Recent Improvements include greater seat depth adjustment (41 cm to 46 cm), a tension adjustable backrest which is extendable from 41 cm to 46 cm and a more robust, rigid castor attachment (0, 3 and 6 degrees).



Scan QR code to
view product online

VARIATIONS

Image	SKU	Price	Description	Size
	S2358	\$799.99		14"
	S9424	\$699.99		16"
	S2351	\$699.99		18"
	S9306	\$699.99		20"

PRODUCT DESCRIPTION

Available in Self Propelled models, the Breezy BasiX2 is a lightweight wheelchair that features amazing adjustability and rigidity, combined with a modern and clean look.

- Adapted seat heights accommodating for foot propulsion
- Armrests conveniently flip back, for easier transfers
- Easy accessible wheel locks prevent the wheelchair from rolling forward or backward
- Angle Adjustable Footplates
- Folding Frame for convenient transportation and storage
- Comes with free standard basic foam cushion
- Swing Away and Adjustable Anti-Tips
- Safe working Load 125 kg
- Fixed Backrest

ADDITIONAL INFORMATION

Brand	Sunrise Medical
Max SWL	125 kg
Seat width	14", 16", 18", 20"
Seat depth	410 mm - 460 mm
Seat height	425 mm - 510 mm
Seat angle	0 degrees to 6 degrees
Back rest height	410 - 460 mm
Seat tilt	No
Frame weight	8.4kg
Warranty	12 months warranty