



**SKU:**

## **SUNRISE MEDICAL BREEZY BASIX 2 FOLDING WHEELCHAIR**





**\$0.00**

The Breezy BasiX<sup>2</sup> is a lightweight folding wheelchair that features amazing adjustability and rigidity, combined with a modern and clean look. Recent Improvements include greater seat depth adjustment (41 cm to 46 cm), a tension adjustable backrest which is extendable from 41 cm to 46 cm and a more robust, rigid castor attachment (0, 3 and 6 degrees).



Scan QR code to  
view product online

VARIATIONS

Image	SKU	Price	Description	Size
	S2358	\$599.99		14 inches
	S9424	\$599.99		16 inches
	S2351	\$599.99		18 inches
	S9306	\$599.99		20 inches

## **PRODUCT DESCRIPTION**

Available in Self Propelled models, the Breezy BasiX2 is a lightweight wheelchair that features amazing adjustability and rigidity, combined with a modern and clean look.

- Adapted seat heights accommodating for foot propulsion
- Armrests conveniently flip back, for easier transfers
- Easy accessible wheel locks prevent the wheelchair from rolling forward or backward
- Angle Adjustable Footplates
- Folding Frame for convenient transportation and storage
- Comes with free standard basic foam cushion
- Swing Away and Adjustable Anti-Tips
- Safe working Load 125 kg
- Fixed Backrest

## ADDITIONAL INFORMATION

<b>Brand</b>	Sunrise Medical
<b>Max SWL</b>	125 kg
<b>Seat width</b>	14", 16", 18", 20"
<b>Seat depth</b>	410 mm - 460 mm
<b>Seat height</b>	425 mm - 510 mm
<b>Seat angle</b>	0 degrees to 6 degrees
<b>Back rest height</b>	410 - 460 mm
<b>Material</b>	Aluminium
<b>Seat tilt</b>	No
<b>Frame weight</b>	8.4kg
<b>Warranty</b>	12 months warranty