



**SKU:**

## **SUNRISE MEDICAL BREEZY BASIX 2 FOLDING WHEELCHAIR**

**\$0.00**

The Breezy BasiX<sup>2</sup> is a lightweight folding wheelchair that features amazing adjustability and rigidity, combined with a modern and clean look. Recent Improvements include greater seat depth adjustment (41 cm to 46 cm), a tension adjustable backrest which is extendable from 41 cm to 46 cm and a more robust, rigid castor attachment (0, 3 and 6 degrees).



Scan QR code to  
view product online

VARIATIONS

| Image                                                                            | SKU   | Price    | Description | Size      |
|----------------------------------------------------------------------------------|-------|----------|-------------|-----------|
|  | S2358 | \$599.99 |             | 14 inches |
|  | S9424 | \$599.99 |             | 16 inches |
|  | S2351 | \$599.99 |             | 18 inches |
|  | S9306 | \$599.99 |             | 20 inches |

## **PRODUCT DESCRIPTION**

Available in Self Propelled models, the Breezy BasiX2 is a lightweight wheelchair that features amazing adjustability and rigidity, combined with a modern and clean look.

- Adapted seat heights accommodating for foot propulsion
- Armrests conveniently flip back, for easier transfers
- Easy accessible wheel locks prevent the wheelchair from rolling forward or backward
- Angle Adjustable Footplates
- Folding Frame for convenient transportation and storage
- Comes with free standard basic foam cushion
- Swing Away and Adjustable Anti-Tips
- Safe working Load 125 kg
- Fixed Backrest

## **ADDITIONAL INFORMATION**

|                         |                        |
|-------------------------|------------------------|
| <b>Brand</b>            | Sunrise Medical        |
| <b>Max SWL</b>          | 125 kg                 |
| <b>Seat width</b>       | 14", 16", 18", 20"     |
| <b>Seat depth</b>       | 410 mm - 460 mm        |
| <b>Seat height</b>      | 425 mm - 510 mm        |
| <b>Seat angle</b>       | 0 degrees to 6 degrees |
| <b>Back rest height</b> | 410 - 460 mm           |
| <b>Material</b>         | Aluminium              |
| <b>Seat tilt</b>        | No                     |
| <b>Frame weight</b>     | 8.4kg                  |
| <b>Warranty</b>         | 12 months warranty     |