




SKU:

## SUNRISE MEDICAL BREEZY EC2000 HD WHEELCHAIR




**~~\$2,499.99~~ \$2,499.99**

A high quality, versatile heavy-duty wheelchair available in seat widths ranging from 22" to 26". The Breezy EC2000 HD wheelchair is a robust, dependable, and versatile product that has your basic mobile needs covered.



Scan QR code to  
view product online

## VARIATIONS

Image	SKU	Price	Description	Size
	S7835	\$999.99		22"
	S7836	\$999.99		24"
	S17165	\$2,499.99		26"

## **PRODUCT DESCRIPTION**

### **Key features**

- Heavy duty reinforced chrome steel frame
- 22"/24" seat widths accommodate users up to 205 kg
- 26" seat widths accommodate users up to 315 kg
- Swing away footrests
- Durable, double-embossed vinyl upholstery (S7835/S7386)
- Durable nylon upholstery (S17165)
- Removable desk length armrests

## ADDITIONAL INFORMATION

<b>Brand</b>	Sunrise Medical
<b>Max SWL</b>	204 kg
<b>Seat width</b>	22", 24"
<b>Seat depth</b>	560 – 760 mm
<b>Backrest height</b>	560 – 760 mm