




SKU: S1306

CUTLERY - CARING SELECTION

\$19.99

The Caring Selection of stainless steel utensils are moulded into contoured, ivory-coloured plastic handles with built-up ends. These are of value to those with a weak grip and restricted wrist or finger movement.

Each sold separately



Scan QR code to
view product online

VARIATIONS

Image	SKU	Price	Description	Type
	S1306	\$19.99		Knife
	S1307	\$19.99		Fork
	S1308	\$19.99		Spoon
	S1309	\$19.99		Angled Spoon RH

PRODUCT DESCRIPTION

The knife and fork have a shaped indent on the top for the index finger, to help with directional control. The knife also has a serrated blade to assist with cutting. In addition to the standard spoon, there is a right handed angled spoon with an angled handle to minimise wrist movement. Hygienically sealed, easy to clean and dishwasher safe.

ADDITIONAL INFORMATION

Brand	ILS ADL
Length of Handle	127 mm
Knife Weight	45 g
Fork Weight	50 g
Spoon Weight	60 g
Angled Spoon Right Weight	60 g
Junior Knife Weight	25 g
Junior Fork Weight	25 g
Junior Spoon Weight	20 g