




SKU: S7267

ROYALE BARIATRIC CHAIR

\$749.99

The Royale Bariatric Highback Chair comes in a 'mocha champagne' colour and has a seating width of 60 cm. The Royale bariatric highback chair has been designed with the user's comfort, posture, and safety in mind. The product features a comfortable, contoured high back design, padded seat and armrests. The chair has height adjustable legs to adjust the seating height for users. The frame has been treated with corrosion resistant powder paint, and the vinyl upholstery is easy to clean and maintain. This product has a safe usage weight of 210 kg and complies with Australian Standards.



Scan QR code to
view product online

PRODUCT DESCRIPTION

The Royale bariatric highback chair is an ideal seating option which provides additional postural support with comfortable padding on the seat, backrest and armrests. The Royale bariatric highback chair comes in a 'mocha champagne' colour and has a seating width of 60 cm. The chair has height adjustable legs to alter seating height, depth adjustable backrest and is designed for multi-purpose use. The chair has durable anti-slip rubber tips, vinyl upholstery and the frame has been treated with corrosion resistant powder paint.

Key Features

- Provides excellent back, upper body, neck, and head support.
- Comfortable padded armrest, seat, and backrest.
- Promotes better posture – easier for user to achieve an upright position.
- Suitable for users requiring assistance to transfer on/off the chair
- Height adjustable seat
- Depth adjustable seat
- Non-slip rubber tips
- Multi-purpose chair, ideally used at the dining table, desk or for general household activities.
- Vinyl upholstery – easy to clean.
- Robust powder coated steel frame.
- Bariatric extra wide seating
- SWL 210 kg

ADDITIONAL INFORMATION

Brand	Royale Medical
Colour	Mocha Champagne
Seat Height (from floor)	460 mm – 610 mm
Seat Depth (measured on sitting area)	480 mm
Seat Width (measured on sitting area)	580 mm
Armrest Height (measured from top of seat)	210 mm
Width (between the arms)	590 mm
Backrest Height (measured from top of seat)	760 mm
Overall Height	1170 mm – 1320 mm
Overall Depth	590 mm
Overall Width	760 mm
Unit Weight	17 kg
SWL	210 kg
ARTG	291700
AS/NZ Standards	AS/NZS 4688.2:2000, 3:2000