




**SKU:**

## **ROYALE MEDICAL HIGH BACK CHAIR**

**\$429.99**

This is the Royale Medical Standard Highback Chair which comes in either a 'mocha champagne', 'forest slate' or 'ardo grey' upholstery. The Royale standard highback chair has been designed with the user's comfort, posture, and safety in mind. The product features a comfortable, contoured high back design, padded seat and armrests. The chair has height adjustable legs to adjust the seating height for users. The frame has been treated with corrosion resistant powder paint and the vinyl upholstery is easy to clean and maintain. This product has a safe usage weight of 205kg and complies with Australian Standards. For additional support, users can purchase the accompanying leg rest (separate purchase).



Scan QR code to  
view product online

### VARIATIONS

Image	SKU	Price	Description	Colour
	S15306	\$429.99		Ardo Grey
	S9705	\$429.99		Mocha Champagne
	S10187	\$429.99		Forest Slate

## **PRODUCT DESCRIPTION**

The Royale Medical standard highback chair is an ideal seating option which provides additional postural support with comfortable padding on the seat, backrest and armrests. The Royale standard highback chair comes in 3 colour variations. The chair has height adjustable legs to alter seating height, depth adjustable backrest and is designed for multi-purpose use. The chair has durable anti-slip rubber tips, vinyl upholstery and the frame has been treated with corrosion resistant powder paint.

### **Key Features**

- Provides excellent back, upper body, neck, and head support.
- Comfortable padded armrest, seat, and backrest.
- Promotes better posture – easier for user to achieve an upright position.
- Suitable for users requiring assistance to transfer on/off the chair
- Height adjustable seat
- Depth adjustable seat
- Non-slip rubber tips
- Multi-purpose chair, ideally used at the dining table, desk or for general household activities.
- Vinyl upholstery – easy to clean.
- Robust powder coated steel frame.

## ADDITIONAL INFORMATION

<b>Brand</b>	Royale Medical
<b>Maximum SWL</b>	205 kg ( 451 lbs)
<b>Seat Height (from floor)</b>	450mm - 600mm
<b>Seat Depth (measured on sitting area)</b>	510mm
<b>Seat Width (measured on sitting area)</b>	510mm
<b>Armrest Height (measured from top of seat)</b>	210mm
<b>Width (between the arms)</b>	490mm
<b>Backrest Height (measured from top of seat)</b>	720mm
<b>Overall Depth</b>	600mm
<b>Overall Width</b>	650mm
<b>Unit Weight</b>	16.4 kg (36.1 lbs)
<b>ARTG</b>	291700
<b>AS/NZS Standards</b>	AS/NZS 4688.2:2000, 3:2000