



**SKU:**

## **ROYALE STANDARD HIGHBACK DAY CHAIR**

**\$429.99**

The Royale Standard High Back Day Chair delivers stable, supportive seating with a contoured backrest and padded design. Adjustable features make it suitable for home, aged care and rehabilitation settings.



Scan QR code to  
view product online

## VARIATIONS

Image	SKU	Price	Description	Colour
	S15306	\$429.99		Ardo Grey
	S9705	\$429.99		Mocha Champagne
	S10187	\$429.99		Forest Slate

## PRODUCT DESCRIPTION

The Royale Standard High Back Day Chair offers ergonomic support in a static seating solution, designed for daily use. Its contoured high backrest promotes optimal posture, while padded armrests assist with standing and transfers.

Seat height and depth can be adjusted to suit individual needs. The chair has anti-slip rubber tips that add safety during use.

The steel frame is powder-coated for durability, and the vinyl upholstery can be cleaned easily.

### Who it's for

- For clients who need reliable, postural support while seated.
- For therapists wanting a comfortable, long-term seating option for their client.

### Key features

- **Contoured high backrest** – supports upper body, neck and head for better posture
- **Padded seat and armrests** – add comfort during long periods of sitting
- **Height and depth adjustable seat** – ensures personalised seating fit
- **Durable vinyl upholstery** – easy to wipe clean and maintain
- **Powder-coated steel frame** – provides strength and resists corrosion
- **Non-slip rubber tips** – enhance safety and stability

## ADDITIONAL INFORMATION

<b>Brand</b>	Royale Medical
<b>Maximum SWL</b>	205 kg ( 451 lbs)
<b>Seat Height (from floor)</b>	450mm - 600mm
<b>Seat Depth (measured on sitting area)</b>	510mm
<b>Seat Width (measured on sitting area)</b>	510mm
<b>Armrest Height (measured from top of seat)</b>	210mm
<b>Width (between the arms)</b>	490mm
<b>Backrest Height (measured from top of seat)</b>	720mm
<b>Overall Depth</b>	600mm
<b>Overall Width</b>	650mm
<b>Unit Weight</b>	16.4 kg (36.1 lbs)
<b>ARTG</b>	291700
<b>AS/NZS Standards</b>	AS/NZS 4688.2:2000, 3:2000