

**SKU:**

# ROYALE STANDARD LOW BACK DAY CHAIR

**\$349.99**



The Royale Standard Low Back Day Chair provides adjustable support and ergonomic seating with a lower backrest. Its contoured design and durable construction make it suitable for various settings, including home and care facilities.



Scan QR code to  
view product online

## VARIATIONS

SKU	Price	Description	Colours
S15305	\$349.99	The Royale Standard Low Back Day Chair provides adjustable support and ergonomic seating with a lower backrest. Its contoured design and durable construction make it suitable for various settings, including home and care facilities.	Ardo Grey
S9706	\$349.99	The Royale Standard Low Back Day Chair provides adjustable support and ergonomic seating with a lower backrest. Its contoured design and durable construction make it suitable for various settings, including home and care facilities.	Forest Slate
S10188	\$349.99	The Royale Standard Low Back Day Chair provides adjustable support and ergonomic seating with a lower backrest. Its contoured design and durable construction make it suitable for various settings, including home and care facilities.	Mocha Champagne

## PRODUCT DESCRIPTION

The Royale Standard Low Back Day Chair offers a supportive and adaptable seating solution for daily use. Its contoured low backrest aids posture and comfort, while padded armrests support ease of transfers.

The chair includes height-adjustable legs and non-slip rubber tips for safer positioning.

Upholstered in durable vinyl, it's easy to clean and maintain across various environments. The steel frame is powder-coated for durability and corrosion resistance.

### Who it's for

- For clients seeking a contoured, Low Back seat for daily use.
- For therapists wanting a practical, low-maintenance chair suited to multiple settings.

### Key features

- **Contoured low backrest** - encourages healthy posture and lumbar support
  - **Padded seat and armrests** - add comfort and assist with transfers
  - **Height-adjustable legs** - ensure a safe, personalised seat height
- **Durable vinyl upholstery** - supports hygiene and easy maintenance
  - **Powder-coated steel frame** - resists wear and corrosion
  - **Non-slip rubber tips** - promote safe and stable use

## **ADDITIONAL INFORMATION**

<b>Safe Working Load (SWL)</b>	205 kg
<b>Overall Width</b>	59 cm
<b>Overall Depth</b>	50 cm
<b>Overall Height</b>	91 cm - 109 cm
<b>Seat Width</b>	46 cm
<b>Seat Depth</b>	46 cm
<b>Seat to Floor Height</b>	47 cm - 65 cm
<b>Backrest Height</b>	45 cm
<b>Brand</b>	Royale Medical

[Royale Standard Low Back User Manual](#)

**Talk to us on the phone**

Our experienced phone sales team can assist you with any product enquiries and help you make a purchase directly on the phone. We are available anytime between 9am and 5pm, Monday to Friday.

[1300 008 267](tel:1300008267)