



EASIGRIP CHIROPODIST SCISSORS

\$60.99

The Chiropodist Scissors long loop handles allow the whole hand to be used for controlling the scissors which is ideal for those with a weak grip. The extended shank helps users who have difficulty bending. The angled blades give a good cutting positions, with the top blade being slightly serrated to stop it slipping. Made from stainless steel blades.



SKU: S1371

Category: [Personal Grooming](#)

ADDITIONAL INFORMATION

Overall Length	200 mm
Blade Length	30 mm
Blade Type	Pointed

Weight

55 g

Brand

ILS ADL