



SKU: S1371

EASIGRIP CHIROPODIST SCISSORS

\$69.99

The Chiropodist Scissors long loop handles allow the whole hand to be used for controlling the scissors which is ideal for those with a weak grip. The extended shank helps users who have difficulty bending. The angled blades give a good cutting positions, with the top blade being slightly serrated to stop it slipping. Made from stainless steel blades.



Scan QR code to
view product online

ADDITIONAL INFORMATION

Overall Length	200 mm
Blade Length	30 mm
Blade Type	Pointed
Weight	55 g
Brand	ILS ADL