

**SKU:** S6954

## COMFORT CHAIR

**\$2,799.99**



The ultimate in comfort and mobility the Comfort Chair is an easy care solution for elderly or special needs seating. Cushioned adjustable air cell support has been incorporated into the design. The four lateral support wings fold away for easy transferring of the patient. This chair has a reclinable seat along with an adjustable backrest and calf rest. The frame is completely constructed from sturdy stainless steel.

Adjustable backrest (90° - 170°) and calf rest

Reclinable seat

Air system for pressure relief

Full stainless steel frame

Three lockable castors and one steering castor (non-marking)

Breathable polyurethane cover (fire retardant)


Replaceable and easy to clean cover

Four folding wings

Standard dining tray

Magazine pocket

10 years warranty on frame



Scan QR code to  
view product online



## **PRODUCT DESCRIPTION**

The comfort chair features an adjustable backrest and calf rest, reclinable seat and air system for pressure relief. A full stainless steel frame with three lockable castors and one steering castor (non-marking) and breathable polyurethane cover (fire retardant) ensures longevity for your comfort air chair. Includes a replaceable and easy to clean cover and four folding wings. Footrest and table come as standard accessories.

## ADDITIONAL INFORMATION

<b>Brand</b>	Days Healthcare
<b>Seat Height</b>	63cm
<b>Seat Width</b>	51cm
<b>Overall Width</b>	73cm
<b>Overall Length</b>	135 - 188cm
<b>Table to Floor</b>	88cm
<b>Table Slide</b>	23cm
<b>Table Size</b>	40 x 40.5 x 2cm
<b>Back Reclining</b>	102° - 136°
<b>Seat Tilting</b>	0 - 18°
<b>Leg Elevating</b>	70° down to 5° up
<b>Weight of Chair</b>	43.5kg
<b>Max SWL</b>	120kg