




SKU:

AIR COMFORT COMPACT LIFT CHAIR

~~\$4,199.99~~ **\$4,199.99**

The KCare Air Comfort Compact Lift Chair is designed to provide pressure relief and support through a specialised air cell system. With lift and recline assistance, it suits individuals with skin integrity concerns or limited mobility. Available in three sizes and with either single or twin motors.



Scan QR code to
view product online

VARIATIONS

SKU	Price	Description	Sizes
S2430	\$3,999.99	The KCare Air Comfort Compact Lift Chair is designed to provide pressure relief and support through a specialised air cell system. With lift and recline assistance, it suits individuals with skin integrity concerns or limited mobility. Available in three sizes and with either single or twin motors.	Single Motor Small
S1469	\$3,999.99	The KCare Air Comfort Compact Lift Chair is designed to provide pressure relief and support through a specialised air cell system. With lift and recline assistance, it suits individuals with skin integrity concerns or limited mobility. Available in three sizes and with either single or twin motors.	Single Motor Medium
S7527	\$3,999.99	The KCare Air Comfort Compact Lift Chair is designed to provide pressure relief and support through a specialised air cell system. With lift and recline assistance, it suits individuals with skin integrity concerns or limited mobility. Available in three sizes and with either single or twin motors.	Single Motor Large

SKU	Price	Description	Sizes
S2236	\$4,199.99	The KCare Air Comfort Compact Lift Chair is designed to provide pressure relief and support through a specialised air cell system. With lift and recline assistance, it suits individuals with skin integrity concerns or limited mobility. Available in three sizes and with either single or twin motors.	Twin Motor Small
S1039	\$4,199.99	The KCare Air Comfort Compact Lift Chair is designed to provide pressure relief and support through a specialised air cell system. With lift and recline assistance, it suits individuals with skin integrity concerns or limited mobility. Available in three sizes and with either single or twin motors.	Twin Motor Medium
S3190	\$4,199.99	The KCare Air Comfort Compact Lift Chair is designed to provide pressure relief and support through a specialised air cell system. With lift and recline assistance, it suits individuals with skin integrity concerns or limited mobility. Available in three sizes and with either single or twin motors.	Twin Motor Large

PRODUCT DESCRIPTION

The KCare Air Comfort Compact Lift Chair is ideal for individuals with pressure care needs, rehabilitation goals or limited mobility. It features the patented Air Comfort seating system, designed to help prevent pressure injuries by redistributing pressure across the seating surface.

The system includes four individual air bladders in the backrest, lumbar, seat and legrest sections, pre-inflated to optimal levels and adjustable if required. Upholstered in VersoTela Fusion material, the chair offers soft, low-shear support suitable for sensitive skin.

With electric lift and recline functions, the Air Comfort Compact makes sitting and standing easier. The seat and backrest cushions are removable for maintenance, while front and rear castors aid with chair positioning.

It's available in three sizes and in single or twin motor options to suit individual needs.

Who it's for

- For clients with pressure sensitivity or skin integrity concerns who need assistive lift and recline support.
- For therapists wanting a clinically designed chair with adjustable air cushions for long-term seating comfort.

Key features

- **Air Comfort seating system** - uses four air cushions to distribute pressure and enhance

comfort

- **Adjustable air bladders** – located in seat, lumbar, backrest and legrest sections
- **VersoTela Fusion upholstery** – low-shear fabric ideal for pressure care and sensitive skin
- **Electric lift and recline functions** – support safe transitions and adjustable positioning
 - **Removable seat and backrest cushions** – simplify cleaning and maintenance
 - **Front and rear castors** – assist with positioning in the home or care setting
- **Choice of single or twin motor** – accommodates individual and clinical preferences

ADDITIONAL INFORMATION

Safe Working Load (SWL)	140 kg
Overall Width	Small: 75 cm, Medium: 75 cm, Large: 82 cm
Overall Depth	Small: 82 cm, Medium: 82 cm, Large: 89 cm
Overall Height	106 cm
Seat Width	Small: 50 cm, Medium: 50 cm, Large: 60 cm
Seat Depth	Small: 44 cm, Medium: 51 cm, Large: 51 cm
Seat to Floor Height	Small: 48 cm, Medium: 52 cm, Large: 52 cm
Number of Motors	Single or Twin
Sizes	Single Motor Small, Single Motor Medium, Single Motor Large, Twin Motor Small, Twin Motor Medium, Twin Motor Large
Brand	KCare Healthcare

[KCare Air Comfort Compact Comparison Guide](#)
[KCare Air Comfort Compact User Manual](#)

Talk to us on the phone

Our experienced phone sales team can assist you with any product enquiries and help you make a purchase directly on the phone. We are available anytime between 9am and 5pm, Monday to Friday.

[1300 008 267](#)

[SINGLE VS TWIN VS QUAD MOTOR LIFT CHAIRS - WHAT'S THE DIFFERENCE?](#)