


SKU:

AIR COMFORT COMPACT LIFT CHAIR

~~\$4,199.99~~ **\$4,199.99**

The KCare Air Comfort Compact Lift Chair is designed to provide pressure relief and support through a specialised air cell system. With lift and recline assistance, it suits individuals with skin integrity concerns or limited mobility. Available in three sizes and with either single or twin motors.



Scan QR code to
view product online



VARIATIONS

Image	SKU	Price	Description	Sizes
	S2430	\$3,999.99	Single Motor Small	
	S1469	\$3,999.99	Single Motor Medium	
	S7527	\$3,999.99	Single Motor Large	
	S2236	\$4,199.99	Twin Motor Small	

Image	SKU	Price	Description	Sizes
	S1039	\$4,199.99	Twin Motor Medium	
	S3190	\$4,199.99	Twin Motor Large	

PRODUCT DESCRIPTION

The KCare Air Comfort Compact Lift Chair is ideal for individuals with pressure care needs, rehabilitation goals or limited mobility. It features the patented Air Comfort seating system, designed to help prevent pressure injuries by redistributing pressure across the seating surface.

The system includes four individual air bladders in the backrest, lumbar, seat and legrest sections, pre-inflated to optimal levels and adjustable if required. Upholstered in VersoTela Fusion material, the chair offers soft, low-shear support suitable for sensitive skin.

With electric lift and recline functions, the Air Comfort Compact makes sitting and standing easier. The seat and backrest cushions are removable for maintenance, while front and rear castors aid with chair positioning.

It's available in three sizes and in single or twin motor options to suit individual needs.

Who it's for

- For clients with pressure sensitivity or skin integrity concerns who need assistive lift and recline support.
- For therapists wanting a clinically designed chair with adjustable air cushions for long-term seating comfort.

Key features

- **Air Comfort seating system** – uses four air cushions to distribute pressure and enhance comfort
- **Adjustable air bladders** – located in seat, lumbar, backrest and legrest sections
- **VersoTela Fusion upholstery** – low-shear fabric ideal for pressure care and sensitive skin
- **Electric lift and recline functions** – support safe transitions and adjustable positioning
- **Removable seat and backrest cushions** – simplify cleaning and maintenance
- **Front and rear castors** – assist with positioning in the home or care setting
- **Choice of single or twin motor** – accommodates individual and clinical preferences

ADDITIONAL INFORMATION

Brand	Air Comfort
Max SWL	140 kg
Back Height	Small: 680 mm / Medium: 740 mm / Large: 800 mm
Seat Height	Small: 480 mm / Medium: 520 mm / Large: 520 mm
Seat Depth	Small: 440 mm / Medium: 510 mm / Large: 510 mm
Seat Width	Small: 500 mm / Medium: 500 mm / Large: 600 mm
Overall Width	Small: 750 mm / Medium: 750 mm / Large: 820mm
Overall Height	1060 mm
Manufacturer's part number	AC59018
Sizes	Single Motor Small, Single Motor Medium, Single Motor Large, Twin Motor Small, Twin Motor Medium, Twin Motor Large
Warranty	A warranty of 5 years on the frame, 2 years on castors and motor 1 year on fabric and struts.