



SKU: S328484

CONFIGURA BARIATRIC 18" HEIGHT BLACK VINYL LIFT CHAIR (TWIN MOTOR)

\$11,499.99

The Configura Bariatric 18" Seat Height Black Vinyl Lift Chair (Twin Motor) supports safe and comfortable seating for individuals with higher body weight. Its 18" seat height, twin-motor function and pressure care features allow personalised adjustment and clinical flexibility.



Scan QR code to
view product online

VARIATIONS

Image	SKU	Price	Description	Size
	S328486	\$11,499.99		24"
	S328487	\$11,499.99		26"
	S328271	\$11,499.99		28"

PRODUCT DESCRIPTION

The Configura Bariatric 18" Seat Height Black Vinyl Lift Chair (Twin Motor) is a clinically focused lift chair designed to meet bariatric seating needs with Configura's advanced adjustability. With independent motors for the backrest and legrest, it enables tailored positioning for posture and comfort.

Its tilt-in-space functionality aids in pressure redistribution, while the lift mechanism assists users in transitioning between sitting and standing positions safely.

The chair includes pressure care support surfaces, a vapour permeable backrest and seat, and a self-supply cushion overlay for clinical versatility. The black vinyl finish offers a low-maintenance, easy-clean solution.

Who it's for

- For clients who need bariatric support, pressure care and an 18" seated height.
- For therapists wanting a configurable, bariatric chair with safe lift assistance and positioning options.

Key features

- **Twin motor design** - enables independent backrest and legrest adjustment
- **Tilt-in-space function** - promotes pressure redistribution and comfort
- **Seat height 18"** - accommodates individuals needing a higher seat profile
- **High weight capacity (up to 254 kg)** - supports bariatric care requirements
- **No shear recline** - preserves pelvic stability during movement
- **Vapour permeable backrest and seat** - protects skin and improves airflow
- **Self-supply cushion overlay** - accommodates personal seating preferences
- **Interchangeable backrest options** - with adjustable Velcro pillows
- **Black vinyl upholstery** - durable and easy to maintain
- **Integrated depth adjustment** - two setting depth adjustment
- **Castors for mobility** - aids relocation in clinical environments
- **Easy-to-use handset** - allows intuitive position control
- **Battery backup** - maintains operation during power outages
- **Range of accessories available** - enhances positioning and clinical fit

ADDITIONAL INFORMATION

Seat Height	40.6 cm (16")
Seat Width	66 cm (26")
Seat Depth	50.8 – 55.8 cm (20" or 22")
Weight Capacity	254 kg
Leg-rest Weight Capacity	127 kg
Brand	Accora
Size	24", 26", 28"