



SKU:

CONFIGURA COMFORT BLACK EDITION TWIN MOTOR LIFT CHAIR

\$4,999.99

The Configura Comfort Black Vinyl (Twin Motor) is a highly adjustable lift chair that promotes comfort, safety and posture support. Its twin-motor system and pressure care features suit a wide range of mobility and clinical needs.



Scan QR code to
view product online

VARIATIONS

Image	SKU	Price	Description	Size
	S312106	\$4,999.99		Small
	S307549	\$4,999.99		Medium
	S310779	\$4,999.99		Large

PRODUCT DESCRIPTION

The Configura Comfort Black Vinyl (Twin Motor) is a versatile electric lift chair designed for those with complex seating needs. With a twin-motor design, it allows the backrest and legrest to be adjusted independently, making it easier to find the right position for comfort and posture.

Its tilt-in-space functionality aids in pressure redistribution, while the lift mechanism assists users in transitioning between sitting and standing positions safely. The chair's vapour permeable seat and backrest work alongside pressure care seating options to support skin integrity.

A self-supply overlay allows the individual's existing cushion to be fitted with ease. The depth can be adjusted, and multiple positioning accessories are available for individualised support. The easy-to-use handset and castors ensure simple operation and repositioning.

The black vinyl upholstery is sleek and easy to maintain, making the chair suitable for a range of care environments.

Who it's for

- For clients who need adjustable posture support, pressure relief and assistance with transfers.
- For therapists wanting a configurable lift chair with clinical functionality and positioning flexibility for their client.

Key features

- **Twin motor design** – enables independent backrest and legrest adjustment
- **Tilt-in-space function** – supports pressure redistribution and comfort
- **No shear recline** – helps maintain pelvic alignment during backrest movement
- **Interchangeable backrest options** – includes adjustable Velcro pillows for customisable support
- **Black vinyl upholstery** – sleek, durable and easy to maintain
- **Vapour permeable seat and backrest** – promotes airflow and protects skin integrity
- **Self-supply cushion overlay** – simplifies fitting of personal seating options
- **Integrated depth adjustment** – two setting depth adjustment
- **Easy-to-use handset** – simplifies position changes
- **Castors for mobility** – enables smooth movement and re-positioning
- **Battery backup** – ensures operation during power outages
- **Available accessories** – supports a range of positioning and clinical needs

ADDITIONAL INFORMATION

Motor	Dual Motor, Twin Motor
Brand	Accora