



SKU:

CONFIGURA TWIN MOTOR WITH TILT RECLINER LIFT CHAIR - BEIGE VINYL

\$4,799.99

The Configura Comfort chair is designed to provide comfortable seating to users with mild to moderate postural requirements. Providing an unprecedented level of functionality and cost-effectiveness.



Scan QR code to
view product online

VARIATIONS

Image	SKU	Price	Description	Size
	S15442	\$4,799.99		Small
	S15443	\$4,799.99		Medium
	S15444	\$4,799.99		Large

PRODUCT DESCRIPTION

Key Features

- Beige Vinyl Upholstery
- Twin Motor
- Tilt-in-space function
- Interchangeable backrest options – adjustable pillow backrests
- Lateral supports and adjustable lateral support
- Vapour Permeable Seat and Backrest
- Easy to use handset
- Pressure Care Seating Options
- Fully hoist accessible with Return Platform and most Stand Hoists, to facilitate reduced carer packages
- Integrated self-supply overlay to make fitting your own cushion easy and simple
- No-shear backrest recline to retain pelvic stability
- Castors for easy mobility
- Knock-down design for easy deliveries
- Integrated depth adjustment – no extra parts required

The Configura Comfort chair is appropriate for:

- Users at risk of pressure damage
- Users with current pressure damage (subject to professional assessment)
- Ambulant, semi-ambulant and non-ambulant users
- Accommodation and encouragement of good posture
- Users with neurological conditions
- Users with a Kyphosis
- Users with oedematous legs and/or fluid build-up

Single vs twin vs quad motor lift chairs - what's the difference?

All ILS lift chairs have a lifting action, to support the user to stand from a seated position.

- Single motor chairs move the back rest and leg rest at the same time.
- Twin motor chairs allow the backrest and leg rest to move independently of each other, typically allowing for a greater degree of leg elevation and greater degree of recline compared to a

standard single motor chair.

- Quad motor chairs allow independent movement of headrest, lumbar, backrest and footrest allowing for the greatest level of re-positioning of all chairs.

ADDITIONAL INFORMATION

Brand	Configura
Max SWL	Small - 130Kg Medium - 160kg Large - 160kg
User Weight Limit	160 kg (352 lbs, 25 Stones)
Seat height	457 mm (18") 406 mm (16")
Seat width	508 mm (20") 508 mm (20")
Armrest height	152 mm (6") 152 mm (6")
Overall width	762 mm (30") 762 mm (30")
Seat length (reclined)	1702mm
Size	Small, Medium, Large