



**SKU:**

## **CONFIGURA TWIN MOTOR WITH TILT RECLINER LIFT CHAIR**




**\$4,499.99**

The Configura Comfort chair is designed to provide comfortable seating to users with mild to moderate postural requirements. Providing an unprecedented level of functionality and cost-effectiveness.



Scan QR code to  
view product online

## VARIATIONS

Image	SKU	Price	Description	Size
	S15442	\$4,499.99		Small
	S15443	\$4,499.99		Medium
	S15444	\$4,499.99		Large

## PRODUCT DESCRIPTION

### Key Features

- Beige Vinyl Upholstery
- Twin Motor
- Tilt-in-space function
- Interchangeable backrest options – adjustable pillow backrests
- Lateral supports and adjustable lateral support
- Vapour Permeable Seat and Backrest
- Easy to use handset
- Pressure Care Seating Options
- Fully hoist accessible with Return Platform and most Stand Hoists, to facilitate reduced carer packages
- Integrated self-supply overlay to make fitting your own cushion easy and simple
- No-shear backrest recline to retain pelvic stability
- Castors for easy mobility
- Knock-down design for easy deliveries
- Integrated depth adjustment – no extra parts required

The Configura Comfort chair is appropriate for:

- Users at risk of pressure damage
- Users with current pressure damage (subject to professional assessment)
- Ambulant, semi-ambulant and non-ambulant users
- Accommodation and encouragement of good posture
- Users with neurological conditions
- Users with a Kyphosis
- Users with oedematous legs and/or fluid build-up

### Single vs twin vs quad motor lift chairs - what's the difference?

All ILS lift chairs have a lifting action, to support the user to stand from a seated position.

- Single motor chairs move the back rest and leg rest at the same time.
- Twin motor chairs allow the backrest and leg rest to move independently of each other, typically allowing for a greater degree of leg elevation and greater degree of recline compared to a

standard single motor chair.

- Quad motor chairs allow independent movement of headrest, lumbar, backrest and footrest allowing for the greatest level of re-positioning of all chairs.

## ADDITIONAL INFORMATION

<b>Max SWL</b>	Small - 130Kg Medium - 160kg Large - 160kg
<b>User Weight Limit</b>	160 kg (352 lbs, 25 Stones)
<b>Seat height</b>	457 mm (18") 406 mm (16")
<b>Seat width</b>	508 mm (20") 508 mm (20")
<b>Armrest height</b>	152 mm (6") 152 mm (6")
<b>Overall width</b>	762 mm (30") 762 mm (30")
<b>Seat length (reclined)</b>	1702mm
<b>Size</b>	Small, Medium, Large