




SKU: S165

BATH BOARD - COSBY

\$299.99

The Cosby bath board provides a stable seating area for bathing and showering. A contoured channeled surface drains water away from seating area while also being reinforced with marine ply timber. Four independently adjustable feet allow the Cosby bathboard to suit most baths. Moulded plastic construction ensures comfort and longevity.



Scan QR code to
view product online

PRODUCT DESCRIPTION

This sturdy bathing assist aids folks who have trouble getting in and out of the tub.

The Cosby Bathboard is divided into two sections. The top is a soft plastic seating area that is free of slats or holes for added comfort. This is reinforced underneath with waterproof plywood, resulting in a bathboard that is both extremely sturdy and extremely slender, measuring only 1 inch in height (25mm).

ADDITIONAL INFORMATION

Width	710 mm
Depth	290 mm
Max SWL	180 kg
daily-living	Bath Aids
Brand	Hero Medical