




**SKU:** S165

## **BATH BOARD - COSBY**

**\$299.99**

The Cosby bath board provides a stable seating area for bathing and showering. A contoured channeled surface drains water away from seating area while also being reinforced with marine ply timber. Four independently adjustable feet allow the Cosby bathboard to suit most baths. Moulded plastic construction ensures comfort and longevity.



Scan QR code to  
view product online

## **PRODUCT DESCRIPTION**

This sturdy bathing assist aids folks who have trouble getting in and out of the tub.

The Cosby Bathboard is divided into two sections. The top is a soft plastic seating area that is free of slats or holes for added comfort. This is reinforced underneath with waterproof plywood, resulting in a bathboard that is both extremely sturdy and extremely slender, measuring only 1 inch in height (25mm).

## ADDITIONAL INFORMATION

<b>Width</b>	710 mm
<b>Depth</b>	290 mm
<b>Max SWL</b>	180 kg
<b>daily-living</b>	Bath Aids
<b>Brand</b>	Hero Medical