

SKU: S12938

COBALT CRASH MAT

\$299.99

The Cobalt Crash Mat provides floor-level protection for individuals at risk of falling from bed. With a high-density foam core and non-slip base, it supports safety in clinical and home care environments.



Scan QR code to
view product online



PRODUCT DESCRIPTION

The Cobalt Crash Mat is a safety solution designed to reduce the impact of accidental bed falls. Positioned beside low or floorline beds, it helps protect individuals who may roll or slip from bed, offering peace of mind in high-risk care scenarios.

The mat is made from high-density foam to absorb impact forces and is finished with a non-slip underside to ensure stability. Its wipe-clean vinyl surface supports infection control, and the foldable design allows easy storage and relocation between care areas.

Who it's for

- For clients at risk of bed related falls who need an added layer of safety.
- For therapists managing clients with dementia, balance concerns or nocturnal movement who require falls prevention strategies.

Key features

- **High-density foam** – cushions falls to reduce injury risk
- **Non-slip base** – maintains position beside the bed
- **Foldable design** – convenient for storage and transport
- **Wipe-clean vinyl cover** – easy to maintain and disinfect
- **Compatible with low beds** – suitable for use with profiling and floorline beds

ADDITIONAL INFORMATION

Brand	Cobalt Healthcare
Thickness	50 mm
Width	600 mm
Length	1750 mm
ARTG	176837
Warranty	2 Years
Folded Dimension	600 x 600 x 150 mm
Mattress Type	Static Foam
Powered or Non-Powered	Non-Powered