

SKU: S12938

CRASH MAT SAFETY MAT

\$299.99



The BetterLiving Safety Crash Mat can be placed on either side of the bed to help prevent injury for patients who roll out of bed. Ideally suited for use with beds in a very low setting, and especially useful when bed safety rails are not suitable.

The unique non-slip feature helps to minimise the risk of injury should the mat be inadvertently stood on by patient or carer. Unique textured underside prevents the mat slipping when downward pressure is applied however still allows the mat to be repositioned with the foot for easy manoeuvrability.


Comprises three sections of safety foam covered by hinged nylon for easy folding and storage, and includes a reflective strip to minimise trip hazards in low light visibility.

helps to prevent patient injury due to falls from bed

sectioned and hinged for easy fold and storage

includes slip resistant textured underside to minimise slip hazard

reflective strip for easy night-time visibility



Scan QR code to
view product online

PRODUCT DESCRIPTION

This Safety Crash Mat can be placed on either side of the bed as a protective cushion to help prevent injury for patients who roll out of bed. Especially useful where side rails cannot be used.

The new non-slip feature helps to minimise the risk of injury should the mat be inadvertently stood on by patient or carer. The unique textured underside prevents the mat slipping when downward pressure is applied however still allows the mat to be repositioned with the foot for easy manoeuvrability.

The Crash Mat includes a reflective strip to minimise trip hazards in low light visibility. Comprises three sections of safety foam covered by hinged nylon for easy folding and storage.

ADDITIONAL INFORMATION

Thickness	50 mm
Width	600 mm
Length	1750 mm
ARTG	176837
Warranty	2 Years
Folded Dimension	600 x 600 x 150 mm