

**SKU:** S18645


## **CURA CORDLESS ALARM BED PAD (NO TRANSMITTER)**

**\$249.99**

Guaranteed 1 year life expectancy  
Use with 2719 Cordless Pad Transmitter  
Easy to clean non-allergenic soft PVC  
Vent holes in the neck reduce the risk of  
moisture ingress

### **Features:**

- Simple operation and maintenance
- Affordable price
- Easy to clean
- For use with cura1 2719 Cordless Pad Transmitter
- Virtually no false alarms over the lifetime of the pad
- All bed/chair pads made from non-allergenic soft PVC
- Vent holes in the neck reduce the risk of moisture ingress



Scan QR code to  
view product online



## **ADDITIONAL INFORMATION**

**Dimensions**

760 x 510mm

**Expectancy**

1 Year life expectancy

**Weight Range**

Recommended weight range for reliable operation: 30kg to 200kg