

SKU: S18645

CURA CORDLESS ALARM BED PAD (NO TRANSMITTER)

\$249.99

Guaranteed 1 year life expectancy

Use with 2719 Cordless Pad Transmitter


Easy to clean non-allergenic soft PVC

Vent holes in the neck reduce the risk of
moisture ingress



Features:

- Simple operation and maintenance
- Affordable price
- Easy to clean
- For use with cura1 2719 Cordless Pad Transmitter
- Virtually no false alarms over the lifetime of the pad
- All bed/chair pads made from non-allergenic soft PVC
- Vent holes in the neck reduce the risk of moisture ingress



Scan QR code to
view product online

ADDITIONAL INFORMATION

Dimensions

760 x 510mm

Expectancy

1 Year life expectancy

Weight Range

Recommended weight range for reliable operation: 30kg to 200kg